

Intended Use

This reagent is intended for *in vitro* quantitative determination of Zinc in serum, plasma & urine.

- Nitro PAPS methodology

Clinical Significance

Zinc is needed for the body's defensive (immune) system to properly work. It plays a role in cell division, cell growth, wound healing, and the breakdown of carbohydrates. Zinc is also needed for the senses of smell and taste. Zinc supplements in large amounts may cause diarrhea, abdominal cramps, and vomiting. Decreased levels are found in cirrhosis, lung carcinomas, sickle cell anemia, acute myocardial infarction, renal failure, corticosteroid and oral contraceptive therapy.

Although zinc is an essential requirement for good health, excess zinc can be harmful.

Principle

Zinc in an alkaline medium reacts with Nitro-PAPS to form a coloured complex. Intensity of the colour is directly proportional to the amount of Zinc present in the sample.

Kit Components

| Reagent/Component | Product Code 12011068 | Description |
|-------------------|-----------------------|---|
| Zinc R1 | 1 x 92 mL | Borate Buffer (pH 8.2) 0.3 M Salicylaldehyde 12.5 mM Dimethylglyoxime 1.25 mM |
| Zinc R2 | 1 x 26 mL | Nitro PAPS 0.4 mM Preservatives |

Risk & Safety

Material Safety data sheets (MSDS) will be provided on request

Reagent Preparation

Zinc R1 & R2 Reagents are ready to use.

Reagent Storage and Stability

The sealed reagents are stable upto the expiry date stated on the label, when stored at 2-8°C.

Open Vial Stability

Once opened the reagents are stable up to 30 days if contamination is avoided.

On-board Calibration Stability

Calibration is stable for 7 days.

Reagent Deterioration

Turbidity or precipitation in any kit component indicates deterioration and the component must be discarded. Values outside the recommended acceptable range for the control may also be an indication of reagent instability and associated results are invalid. Sample should be retested using fresh vial of reagent.

Precaution

To avoid contamination, use clean laboratory wares. Close reagent bottles immediately after use. Do not blow into the reagent bottles.

This reagent is only for IVD use and follow the normal precautions required for handling all laboratory reagents.

Waste Management

Reagents must be disposed off in accordance with local regulations.

Sample

Serum/Plasma (Do not use lipemic or hemolysed sample). Use Heparin as anticoagulant.

Interferences

No interference for

Bilirubin up to 20 mg/dL

Materials provided

Zinc R1 & R2 Reagent

Reagents required but not provided

Zinc Calibrator (Product Code: 11640001)

Unit Conversion

| Traditional Unit | SI Unit | Conversion from Traditional to SI |
|------------------|---------|-----------------------------------|
| µg/dL | µmol/L | x 0.153 |

Calibration

Zinc Calibrator (Product Code: 11640001) is recommended for calibration of the assay.

Quality Control

It is recommended to use Control to verify the performance of the measurement procedure.

Each Laboratory has to establish its own internal quality control scheme and procedures for corrective action if controls do not recover within the acceptable tolerance.

Reference Range

It is recommended that each laboratory should establish its own reference values. The following value may be used as guide line.

Serum/Plasma : 70 - 115 µg/dL

Urine (24 hours) : 100 - 1000 µg/dL

Results obtained for patient samples are to be correlated with clinical findings of patient for interpretation and diagnosis.

Performance

1. Linearity

The reagent is linear up to 1000 µg/dL. If the concentration is greater than linearity (1000 µg/dL), dilute the sample with normal saline and repeat the assay. Multiply the result with dilution factor.

2. Comparison

A comparison study has been performed between Agappe reagent and another internationally available reagent yielded a correlation coefficient of $r^2 = 0.998$ and a regression equation of $y = 0.985x$.

3. Precision

| | Intra Run | | Inter Run | |
|--------------|-----------|---------|-----------|---------|
| Control | Level 2 | Level 3 | Level 2 | Level 3 |
| n | 20 | 20 | 20 | 20 |
| Mean (µg/dL) | 147.97 | 215.41 | 142.81 | 219.16 |
| SD | 2.4 | 3.9 | 4.1 | 1.7 |
| CV(%) | 1.6 | 1.8 | 2.9 | 0.8 |

Accuracy (µg/dL)

| Control | Expected Value | Measured Value |
|-----------------|----------------|----------------|
| Control Level 2 | 127 ± 25 | 139.0 |
| Control Level 3 | 219 ± 43 | 219.0 |

4. Sensitivity

Lower detection Limit is 29 µg/dL

Bibliography

- Pasquinelli, F.; Diagnostica e Tecniche di Laboratorio, (pag.:1103-1104) Rossini Editrice.(1984)
- Testuuo Makino; Chimica Clinica Acta 197, 209-220(1991)
- Maringonia, A., Illuzi, R, ATB 1991 Abstract.

SYMBOLS USED ON THE LABELS

IVD IN VITRO DIAGNOSTIC USE SEE PACKAGE INSERT FOR PROCEDURE LOT LOT NUMBER MANUFACTURER'S ADDRESS MANUFACTURING DATE EXPIRY DATE TEMPERATURE LIMIT



AGAPPE DIAGNOSTICS LTD.

'Agappe Hills', Dist. Ernakulam, Kerala, India-683 562.
Tel. +91 484 2867 000 | Customer Support No.: 1800 425 7151 (Toll Free)
customersupport@agappe.in | www.agappe.com

REV. NO.: ADL/IFU/ZN/120FR/R02

CE ISO 9001:2015
EN ISO 13485:2016