

## CALCIUM (ARSENAZO)

2 x 30 mL  
12013007

### Intended Use

This reagent is intended for *in vitro* quantitative determination of Calcium in serum, plasma & urine.

- Modified Arsenazo III method
- Linear up to 16 mg/dL

### Clinical Significance

Calcium is an important ion present in the body. It is mainly found in bones. In serum calcium exists equally in a free ionized form & also in a bound form with albumin. Calcium helps in enzyme activation, muscle contraction, coagulation of blood, regulation of some hormonal secretions & cell membrane permeability. Increased levels are found in hyperthyroidism, malignant tumors, acute & osteoporosis, adrenal insufficiency. Decreased levels are found in hypoparathyroidism, osteomalacia, rickets, renal failure & tetanus.

### Principle

At a neutral pH the  $Ca^{2+}$  form with Arsenazo III a complex, the color intensity of which is directly proportional to the concentration of calcium in the sample.

### Kit Components

Reagent/Component	Product Code	Description
Calcium Arsenazo Reagent	12013007	MES, pH 6.5 Arsenazo III
		1000 mmol/L 200 mmol/L

### Risk & Safety

Material Safety data sheets (MSDS) will be provided on request

### Reagent Preparation

Calcium Arsenazo reagent is ready to use.

### Reagent Storage and Stability

The sealed reagents are stable upto the expiry date stated on the label, when stored at 2-8°C.

### Open Vial Stability

Once opened the reagents are stable up to 14 days when contaminations are avoided.

### Reagent Deterioration

Turbidity or precipitation in any kit component indicates deterioration and the component must be discarded. Values outside the recommended acceptable range for the Agappe Qualicheck Norm & Path control may also be an indication of reagent instability and associated results are invalid. Sample should be retested using fresh vial of reagent.

### Precaution

To avoid contamination, Close reagent bottles immediately after use. Avoid direct exposure of reagent to light. Do not blow into the reagent bottles. This reagent is only for IVD use and follow the normal precautions required for handling all laboratory reagents.

### Waste Management

Reagents must be disposed off in accordance with local regulations.

### Sample

Fresh serum / plasma (Do not use lipemic or hemolysed sample)  
Urine diluted 1/3 with distilled water; adjust to pH 3-4 with HCl (N/10).  
Take dilution factor into account for the calculation of the concentration in urine.

### Interferences

No interference for  
Bilirubin up to 10 mg/dL  
Haemoglobin up to 1000 mg/dL  
Ascorbic acid 50 mg/dL

### Materials provided

Calcium Arsenazo reagent

### Reagents required but not provided

Multicalibrator (Product Code: 11610001), Qualicheck Norm (Product Code: 11601003), Qualicheck Path (Product Code: 11601002)

### Unit Conversion

Traditional Unit	SI Unit	Conversion from Traditional to SI
mg/dL	mmol/L	x 0.25

### Calibration

Agappe Multicalibrator (Product Code: 11610001) is recommended for calibration of the assay.

### Quality Control

It is recommended to use Qualicheck Norm (Product Code: 11601003) or Qualicheck Path (Product Code: 11601002) to verify the performance of the measurement procedure. Each Laboratory has to establish its own internal quality control scheme and procedures for corrective action if controls do not recover within the acceptable tolerance.

### Reference Range

It is recommended that each laboratory should establish its own reference values. The following value may be used as guide line.

Serum/plasma: 8.8 - 10.2 mg/dL

Urine: 100 - 400 mg/24 hrs

Results obtained for patient samples are to be correlated with clinical findings of patient for interpretation and diagnosis.

### Performance

#### 1. Linearity

The reagent is linear up to 16 mg/dL. If the concentration is greater than linearity (16 mg/dL), dilute the sample with normal saline and repeat the assay. Multiply the result with dilution factor.

#### 2. Comparison

A comparison study has been performed between Agappe reagent and another internationally available reagent yielded a correlation coefficient of  $r^2 = 0.9791$  and a regression equation of  $y = 0.9951x$ .

#### 3. Precision

	Intra Run		Inter Run	
Control	Level 1	Level 2	Level 1	Level 2
n	20	20	20	20
Mean (mg/dL)	9.8	12.6	9.81	12.71
SD	0.18	0.15	0.16	0.17
CV(%)	1.85	1.22	1.60	1.31

Accuracy (mg/dL)		
Control	Expected Value	Measured Value
Control Level 1	9.47 ± 0.93	9.7
Control Level 2	12.5 ± 1.3	12.2
Qualicheck Norm	8.80 ± 0.82	9.0

### 4. Sensitivity

Lower detection Limit is 0.5 mg/dL

### Bibliography

1. Bayer, P.J. Anal. Biochem., 110, (1981), 61

### Mispa Mini Assay Parameter

Page 1	
Test (Name)	Calcium A
Decimal	2
Unit	mg/dL
Reac. Method	End point
Blank	R blank
Wave length I	630
Wave length II	-
Test Time	300 - 300
Sample	3 µL
R1	300 µL
Calculation Method	One Point Linear

### SYMBOLS USED ON THE LABELS

IVD IN VITRO DIAGNOSTIC USE SEE PACKAGE INSERT FOR PROCEDURE LOT LOT NUMBER MANUFACTURER'S ADDRESS MANUFACTURING DATE EXPIRY DATE TEMPERATURE LIMIT



**AGAPPE DIAGNOSTICS LTD.**

'Agappe Hills', Dist. Ernakulam, Kerala, India-683 562.  
Tel. +91 484 2867 000 | Customer Support No.: 1800 425 7151 (Toll Free)  
customersupport@agappe.in | www.agappe.com

REV. NO.: ADL/IFU/CALA/MINI/R01



ISO 9001:2015  
EN ISO 13485:2016